
LASOS[®] DPSS laser series

*LK wwx Tyzz

Diodengepumpter Festkörperlaser
Diode Pumped Solid-State Laser

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It is a requirement to know the contents of this manual before starting up the laser system. Therefore, please read this manual carefully and pay special attention to the notes regarding the safe use of the system.

This manual describes the use of the standard system for stand alone use. The laser system is also available as OEM versions with different configurations. In this case the system is clearly marked as OEM component and the operation may differ.

© LASOS and the LASOS logo are registered trademarks of LASOS Lasertechnik GmbH.

This document is copyrighted with all rights reserved. This also applies in the event of granting of patents and registration of utility models. This document may not be copied in whole or in part unless LASOS has given its express written permission. This includes reproduction or publication in any other media as well as the translation into other languages.

These operating instructions have been made with great care to ensure that data and information given in this document are correct. LASOS makes no warranty or representation neither expressed nor implied with respect to this document. Especially, LASOS will in no event be liable for any direct or indirect damage resulting from defects in this document.

Modifications of drawings, tables, specifications and other information are subject to change without prior notice. These operating instructions are not covered by an update service.

Further information can be obtained direct from the manufacturer or your national distributor. You will find the up-to-date list of our distributors on the Internet under <http://www.lasos.com>.

Manufacturer's address:	LASOS Lasertechnik GmbH Franz-Loewen-Str. 2 07745 Jena Germany
Phone:	(+49) 3641 / 29 44-0
Fax:	(+49) 3641 / 29 44-300
Internet:	http://www.lasos.com info@lasos.com

Number of the operating instructions: 577039-5133-000

Date of publication: 0B / 20.August 2013

Contents

	Page
1	Safety instructions.....5
1.1	General / Symbols 5
1.2	Electrical safety / Protection against accidental contact 6
1.3	Laser safety 7
1.4	Waste removal 8
1.5	Notes on warranty 9
1.6	Intended use 9
2	Device name / device models10
2.1	Overview about allowed Stand-Alone configurations..... 10
2.2	Equipment with controller CKL 10xy z 11
2.3	Equipment with controller CKL 24xy z 12
3	Operation.....13
3.1	Starting up 13
3.1.1	Unpacking / installation 13
3.1.2	Configuration with CKL 24xy z..... 15
3.1.2.1	Controller CKL 24xy z 16
3.1.3	Configuration with CKL 10xy z..... 19
3.1.3.1	CKL 10xy z 20
3.1.3.2	SwitchBox SSL 1000 21
3.1.4	Connection Cable CKL1000..... 23
3.1.5	Power supply 25
3.1.6	Power connection 25
3.2	Switching the system on/off 26
3.2.1	Switching on (CKL 24xy z)..... 26
3.2.2	Switching on (CKL 10xy z)..... 27
3.2.3	Switching off (CKL 24xy z / CKL 10xy z) 29
4	Maintenance30
4.1	Care..... 30
4.2	Servicing 30
4.3	Transportation 31
5	Technical data.....32
5.1	General technical data 32
5.2	Laser types and useable accessories..... 33
5.3	Dimension drawings: 36
5.4	Maximum thermal output of the base plate of the laser module..... 38
6	LASOS Communicator.....39
6.1	General 39
6.2	Getting started..... 39
6.2.1	Installation of the software..... 39

6.2.2	Connection to the laser system.....	40
6.2.2.1	Connecting by USB	40
6.2.2.2	Connecting by RS232.....	40
6.3	Using the LASOS Communicator.....	41
 EC declaration of conformity.....		43

1 Safety instructions

1.1 General / Symbols

Please read these instructions carefully before starting up the device.

For questions or further information, please contact our service team or visit us on the Internet at www.lasos.com.

The LASOS® DPSS laser series is designed, manufactured and tested complying with the standards DIN EN 61010-1 (IEC 61010-1 "Safety requirements for electrical equipment for measurement, control and laboratory use").

The LASOS® DPSS laser series complies with the Directive 2004/108/EC of the European Parliament and of the Council of 15 December 2004.

The LASOS® DPSS laser series is labelled with **CE** - mark.

The LASOS® DPSS laser series will be disposed of according to the WEEE Directive 2012/19/EC of the European Parliament and of the Council of 04 July 2012.

These operating instructions contain the following information and warning symbols:

NOTE

This symbol refers you to notes containing important information on operation.

WARNING

This symbol refers you to a specific procedure or practice which, if not followed correctly, could cause personal injury.

CAUTION

This symbol refers you to a procedure which, if not observed, could result in damage to the equipment.

WARNING

Emission of laser radiation!

CAUTION

Components sensitive to electrostatic discharge.

CAUTION

Take appropriate precautions to prevent damage to components sensitive to electrostatic discharge (packaging, grounding wrist strap, or similar things).

1.2 Electrical safety / Protection against accidental contact

The laser module may only be operated in combination with approved accessory (see the listing of approved accessory on page 33), which is to be connected to a power supply that is galvanically isolated from line power.

The laser system has to arrange that the mains plug can pull out in any situations without problems.

All cables must be sufficiently isolated. Only service technicians employed with or authorized by the manufacturer are allowed to open the devices.

The touch protection of the laser module must be ensured by the user.

To prevent damage to the device, never connect or disconnect the cable between the Controller and the laser module while the power supply is switched on.

1.3 Laser safety

The laser system in combination with permitted Controllers and SwitchBox complies to the standards of CDRH (CFR 21 1040.10) and EN 60825-1.

Laser safety of the laser module must be ensured by the user.

Never open the case of the laser module! Inside the laser module, dangerous laser radiation exists.

For installation and operation, observe the following relevant regulations: DIN EN 61010-1:2011-07, DIN EN 60950:2012-05, DIN EN 60825-1:2012-11 and BGV B2.

After you have switched on the power supply, laser radiation can be emitted by the connected laser module. This state is indicated by lighting of the yellow emission indicator lamp located at the front panels of the SwitchBox and red emission indicator lamp of the laser module.

If provided, the key switch serves to switch the laser radiation on and off when the power supply is switched on. An authorized person should keep the key (Fig. 19/1, Fig. 20/1) in a safe place.

Post a laser warning sign at the place of installation of the laser module so that it is clearly visible.

The safety interlock circuit has to be connected with appropriate safety devices at the interlock connector of the SwitchBox or CKL 2400 xx (if available).

Secure the laser module against unintentional displacement and avoid setup at the user's eye level.


Make sure to avoid uncontrolled reflections by reflecting surfaces or similar objects. In the laser beam area, only use tools with appropriate, i.e. not directly reflecting surfaces.

Before the first start-up of the system, install all protective devices necessary for the safe operation of the laser module. The responsibility for protecting the laser module against improper use by a key switch solely lies with you.

DPSSL

STANDARD :

OEM :


Fig. 1 Warning labels / type labels on system components (examples)

Laser class

Warning! Class IIIb laser device

Avoid exposure of eye or skin to direct or scattered radiation.
Never look directly or by means of optical instruments into the laser beam!

The device complies with the following safety regulations: Center for Devices and Radiological Health (CDRH) Laser Class IIIb, DIN EN 60825-1:2008-05 Laser Class 3B. It is labeled with the corresponding warning signs (see Fig. 1).

The warning notes are part of the system. Always observe these warning notes when you work on it.

Never put the system into operation if warning notes are damaged or missing and immediately inform the competent person (laser safety officer) accordingly.

1.4 Waste removal

Do not dispose of unserviceable equipment with normal household waste. They have to be disposed of according to legal requirements.

The manufacturer of the equipment is legally obligated to take unserviceable equipment back.

1.5 Notes on warranty

The manufacturer warrants that the device has been supplied in a state free from defects of material and/or workmanship. It is imperative that you report any defect to the manufacturer immediately and do everything to minimize the damage.

The warranty does not cover any defects caused by natural wear, especially in case of wearing parts as well as inappropriate use of the device.

The manufacturer is not liable for damages which result from operating errors, negligence or unauthorised tampering with the equipment, particularly as the result of removal or replacement of parts of the device or as the result of the use of unsuitable accessories of other make. Mind that any such action will render all warranty claims null and void.

1.6 Intended use

Do not operate the DPSSL as stand-alone laser without corresponding safety devices.

The DPSSL serve to generate continuous-wave (cw) laser radiation.

The DPSSL operates with output powers ranging from 1 mW to 200 mW.

The operation of the laser module with the approved Controller and SwitchBox (see table at page 33) ensures compliance with the requirements of EC Directive 2004/108/EG and with the Electromagnetic Compatibility Law as amended on February 26, 2008 according to the following standards:

DIN EN 61000-6-2:2006-03

DIN EN 61000-6-3:2011-09

The laser module may only be operated by authorized persons who were instructed in the use of the laser and who know the applicable laser safety regulations.

Keep this manual ready to hand near the device so that it is accessible to all users.
Verify that the necessary documentation, including work instructions with laser-safety regulations or national safety regulations is available, too.

2 Device name / device models

2.1 Overview about allowed Stand-Alone configurations


Fig. 2 Configuration with controller CKL 10xy z



Fig. 3 Configuration with controller CKL 24xy z

2.2 Equipment with controller CKL 10xy z

Manufacturer's designation: LASOS® DPSS laser series

Model name: *LK www Tyzz

- | | |
|---------------------------------|---------------------------------|
| 1 DPSSL | 6 Key for key switch |
| 2 Transport connector (if used) | 7 RS 232 connector |
| 3 Controller CKL 10xy z | 8 Power supply |
| 4 Connection Cable CKL1000 | 9 Power cord (country-specific) |
| 5 SwitchBox SSL 1000 | |

Fig. 4 Overview of system components (with CKL 10xy z)

The components supplied with the laser system may differ at customized versions. Refer to the delivery notes.

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
Fig. 5 Overview of system components (with CKL 24xy z)

- | | |
|---------------------------------|---------------------------------|
| 1 DPSSL | 5 Key for key switch |
| 2 Transport connector (if used) | 6 Power supply |
| 3 Controller CKL 24xy z | 7 Power cord (country-specific) |
| 4 Connection cable CKL 2400 | |

The components supplied with the laser system may differ at customized versions. Refer to the delivery notes.

3 Operation

3.1 Starting up

The DPSSL can be commissioned independently by the user.

Before installing and starting up the system, carefully read the safety instructions at the beginning of this manual.

Before commissioning the laser module, install all safety devices needed for the safe operation of the laser module. The responsibility for protecting the laser module against improper use by a key switch solely lies with you.

3.1.1 Unpacking / installation

The laser diode of the laser module is very sensitive to electrostatic discharge. The discharge of static electricity (e.g. through fingers) can reduce the lifetime of the laser module or cause total damage. To prevent this, take appropriate precautions (grounding wrist strap or a similar aid).

- Please, inspect the shipping boxes for signs of rough handling or damage. Any damage must be reported immediately to shipping carrier, and to manufacturer.
- Remove the laser module and the other units from the packaging and check it for any damage caused by transportation.

- If provided do not remove the transport connector (Fig. 6/1) from the laser-head connector (female connector 24-pin, (Fig. 6/2) unless you connect the control unit.

Usually, the controller device with the Connection Cable Controller is already mounted.

Cooling of the laser module is via its base plate.

- Mount the laser module to a heat sink (e.g. flat aluminum plate, flatness 0.02 mm) ensuring heat dissipation of 12 W.

While in operation, the temperature of the base plate of the laser module must not exceed 45 °C.

For mounting, use the provided bores at the undersides of laser module and control unit. For the mounting dimensions, refer to the Dimension drawing on page 36 or the data-sheet.

Insufficient heat dissipation may cause overheating and total failure of the laser module.

If used the cable is pre-installed between DPSSL and controller.

Only authorized persons using appropriate tools are allowed to connect the laser module to the controller and/or other components like the SwitchBox.



Fig. 6

Remove the transport connector

3.1.2 Configuration with CKL 24xy z

Manufacturer's designation: LASOS® DPSS laser series
Model name laser module: *LK ww x Tyzz
Controller: CKL 24xy z

- | | |
|---------------------------------|---|
| 1 Laser module *LK ww x Tyzz | 4 Key for key switch |
| 2 Transport connector, optional | 5 Laser connection cable |
| 3 Controller CKL 24xy z | 6 Hexagon socket head screw ISO 4762 M3x6, optional |

Fig. 7 Overview of the components of CKL 24xy z

The control unit is available as OEM model in different modifications.

3.1.2.1 Controller CKL 24xy z

The controller allows you to communicate with the laser module via RS232. Please inform about the facilities of connection and communication in the chapter 6 (LASOS Communicator).

Furthermore the controller has the Remote Interlock.

Only authorized persons using appropriate tools are allowed to connect the laser module to the control unit only.

- Connect the connecting cable with its 24-pin female connector (Fig. 8/3) to the male connector on the laser module (Fig. 8/2).
- Connect the other end of the connecting cable with its 15-pin SUB-D male connector (Fig. 8/4) to the Sub-D female connector **LASER** of the control unit (Fig. 8/1).
- Secure both ends of the connecting cable with the two fastening screws against accidental slip off.


Fig. 8 Connecting the laser module

- Connect the power supply to the control unit (Fig. 9/1) by means of male connector KPP-3P (KYCON) (Fig. 9/3). Pay attention to the connected loads from the data sheet of the control device. Connect every laser module used to a separate power supply, if possible.
- Connect the safety interlock circuit to the control unit (if INTERLOCK provided) by means of male connector KMDLA-3P (KYCON) (Fig. 9/2).

The responsibility for the use of a safety interlock circuit lies with the user.


Fig. 9 CKL 2400, connector for operating voltage and interlock

Interlock specifications	Value	
Disconnecting time	<20ms	
Signal status	open close	laser off laser on
Impedance	>1000Ω	
Current Interlock circuit	<10mA	


If the safety interlock circuit connected to the control unit is not closed, the laser does not work. If the safety interlock circuit is interrupted during laser operation, the laser instantly switches off. To restart the laser system, it is necessary to disconnect the power supply by about 1 min before reconnecting it again.

The laser is switched on (Fig. 19/**on**) and off (Fig. 19/**off**) by the key switch (Fig. 19/1).

The emission control lamp (Fig. 19/2) signalises the presence of dangerous laser radiation.

3.1.3 Configuration with CKL 10xy z

Manufacturer's designation: LASOS® DPSS laser series
Model name laser module: *LK wwx Tyzz
Controller: CKL 10xy z
SwitchBox: SSL 1000

- | | |
|---------------------------------|--|
| 1 DPSSL *LK wwx Tyzz | 5 SwitchBox SSL 1000 |
| 2 Transport connector (if used) | 6 Key for key switch |
| 3 Controller CKL 10xy z | 7 Laser connecting cable |
| 4 Connection cable CKL1000 | 8 Hexagon socket head bolt ISO 4762 M3x6; optional |

Fig. 10 Overview of the components of CKL 10xy z

The control unit is available as OEM model in different modifications.

3.1.3.1 CKL 10xy z


Fig. 11 CKL 10xy z Front side / Back side


Fig. 12 Connecting DPSSL

- If not mounted connect the connecting cable with its 24-pin female connector (Fig. 12/3) to the male connector at the laser module (Fig. 12/2).
- Connect the other end of the connecting cable with its 15-pin SUB-D male connector (Fig. 12/4) to the Sub-D female connector **LASER** of the control unit (Fig. 12/1).
- Secure both ends of the connecting cable with the two fastening screws against accidental slip off.

3.1.3.2 SwitchBox SSL 1000

- Connect the the Connection Cable CKL (9-pin D-Sub male connection Fig. 13/1) with the SwitchBox (Connection Cable CKL / 9-pin D-Sub female connector Fig. 13 /2) (if present).


Fig. 13 Install Connection Cable CKL1000 at the SwitchBox

- Connect the safety loop (Interlock, Fig. 14/2) to an appropriate external safety device. Use a 2 mm laboratory plug for connecting the safety loop (Fig. 14/1).


Fig. 14 Safety loop (Interlock)

The responsibility for the use of a safety interlock circuit lies with the user.

The contacts of the Interlock must always been wired with zero potential.

The analyse electronic for the Interlock is provided in the controller CKL 10xy z.

Interlock specifications	Value	
Disconnecting time	<20ms	
Signal status	open	laser off
	close	laser on
Impedance	> 1000Ω	
Current Interlock circuit	<5mA	

If the safety interlock circuit connected to the control unit is not closed, the laser does not work. If the safety interlock circuit is interrupted during laser operation, the laser instantly switches off. To restart the laser system, it is necessary to disconnect the power supply by about 1min before reconnecting it again.

3.1.4 Connection Cable CKL1000


Fig. 15 Connection cable CKL1000

- | | |
|-------------------------|------------------------|
| 1 Connection controller | 4 RS232 connection |
| 2 Power supply | 5 Connection SwitchBox |
| 3 Mains supply | |

The Connection cable CKL1000 connects all components of the laser systems. Furthermore the cable has a D-Sub male connector for the RS232 connector (Fig. 15/6). The electrical power is provided by the fix connected power supply (Fig. 15/2).


Fig. 16 Install Connection Cable CKL1000 at the SwitchBox

- Connect the Connection Cable CKL 1000 (15 pin HD D-Sub-female connector; Fig. 16/1) with the SwitchBox (15-pin HD D-Sub male connector Fig. 16/2).

3.1.5 Power supply

The power supply is a universal AC/DC converter with far range input 12V / 36W and stabilized output (Fig. Fig. 17/2).

Any unauthorized change of the cables may damage the system. It may lead to loss of compliance to the standards listed in chapter 1.

3.1.6 Power connection

- Connect the IEC 60320-mains input connector of the power supply (Fig. Fig. 17/1) with a mains socket (Fig. Fig. 17/3) by means of a power cord (Fig. Fig. 17/2).


Fig. 17 Mains connection

The mains connector of the power cord can vary depending on your country.

3.2 Switching the system on/off

3.2.1 Switching on (CKL 24xy z)

The laser system (laser module and control unit) is switched on by applying a supply voltage to the connected control unit. The red emission indicator lamp on control unit (Fig. 19/2) light up.

Before, an authorized person must insert the key (Fig. 19/1) in the key switch, if used.

- If provided, switch the key switch (Fig. 19/1) of the control unit in the **on** position.

The red emission indicator lamp at the laser module (Fig. 18/1) light up.

The laser reaches its full power after a warm-up time of maximally 15 minutes.

- Open the laser beam path by opening the shutter (if provided) on the laser module (Fig. 18/2). Do not open the beam path unless the working conditions are safe!


Fig. 18 laser emission indicator lamp and shutter

After applying a supply voltage and switch-on of the control unit (if key switch is provided), instantly laser radiation can be emitted from the connected laser module. This condition is signaled by lighting of the red emission indicator lamps located at the top of laser module (Fig. 18/1) and at the front panels of control unit (Fig. 19/2).


Fig. 19 Switching on the control unit

If the safety interlock circuit connected to the control unit (if INTERLOCK provided) is **not closed**, the laser does not work. If the safety interlock circuit is interrupted during laser operation, the laser instantly switches off. To restart the laser system, it is necessary to disconnect the power supply and then reconnect it again.

With the key switch you are able to turn the Laser ON/OFF. The key switch has highest priority in the system.

The laser system provides several control and start possibilities.

- Control via RS232
The RS232 connector is provided at the front of the controller. To communicate with the controller refer to the description in chapter 6 LASOS communicator).

3.2.2 Switching on (CKL 10xy z)

The laser system is switched on by applying a supply voltage to the SwitchBox. The yellow emission indicator lamp on the SwitchBox (Fig. 20/2) lights up.

- If provided, switch the key switch (Fig. 20/1) of the control unit in the **on** position.

The red emission indicator lamp at the laser module (Fig. 18/1) lights up.

An authorized person must insert the key (Fig. 20/1) in the key switch, if used.


Fig. 20 Key switch


Fig. 21 Emission control and shutter

The laser reaches its full power after a warm-up time of maximally 15 minutes.

- Open the laser beam path by opening the shutter (if provided) on the laser module (Fig. 21/2). Do not open the beam path unless the working conditions are safe!

After applying a supply voltage and switch-on of the SwitchBox (if key switch is provided), instantly laser radiation can be emitted from the connected laser module. This condition is signaled by lighting of the red emission indicator lamp located at the top of the laser module (Fig. 18/1) and the yellow emission indicator lamp at the front panels of SwitchBox (Fig. 20/2).

If the safety interlock circuit connected to the SwitchBox is not closed, the laser does not work. If the safety interlock circuit is interrupted during laser operation, the laser instantly switches off. To restart the laser system, it is necessary to disconnect the power supply by about 1 min before reconnecting it again.

3.2.3 Switching off (CKL 24xy z / CKL 10xy z)

The laser system is switched off by disconnecting the voltage from the connected power supply.

Please consider the switch off order!

- Close the laser beam path with the shutter (if provided) on the laser module (Fig. 21/closed).
- Turn the key of the key switch at the connected SwitchBox or the controller CKL 24xy z in **off** position respectively (Fig. 19/off, Fig. 20/off), if provided.
- Turn off the supply voltage.

The laser module is switched off. The yellow emission indicator lamps on SwitchBox (Fig. 20/2) and the red emission indicator lamp on the laser module (Fig. 21/1) go out.

The key is removable in position "0". (Fig. 20/off)

The key (Fig. 20/1) has to be stored at a safe place.

If provided an authorized person shall remove the key (Fig. 19/1, Fig. 20/1) from the device and store it in a safe place for safety reasons.

4 Maintenance

4.1 Care

Before any care, switch off the laser system (see 3.2.3)!

Please make sure to observe the cleaning instructions provided with equipment of other make used in combination with the laser system.

- To clean the device, never use cleaning powders, paint thinners or solvents, such as benzine or acetone, as these attack the surface of the device case.
- Clean the device, if necessary, with a soft cloth moistened with mild soapsuds or a disinfectant.

4.2 Servicing

Service, maintenance, adjustment and cleaning should be performed exclusively by customer service personnel or by specially authorized persons.

4.3 Transportation

Observe the following instructions for safe transportation. Non-observance of these instructions will void all warranties.

The laser diode of the laser module is very sensitive to electrostatic discharge. The discharge of static electricity (e.g. through fingers) can reduce the lifetime of the laser module or cause total damage. If required to prevent this, take appropriate precautions (grounding wrist strap or a similar aid).

Only use packaging material appropriate for components that are sensitive to electrostatic discharge. Avoid the use of styrofoam or simple air-cushion films. For re-packing system components preferably use the original packaging.

- Before any transportation or storage of the device, close the shutter (if provided) at the laser exit aperture (Fig. 21/closed) to avoid contamination.
- If required store and carry the laser diode module with the connected cables (if present).
- Ship the laser with the connecting cable. If only the laser module has to be shipped plug the transport connector (Fig. 22/1) into the laser connector on the laser module (Fig. 22/2).
- Transport the device preferably in its original packaging.
- For the transportation of the laser system observe the ambient conditions specified in chapter 6 Technical data on page 32.


Fig. 22 **Preparing the laser module for transportation / storage**

5 Technical data

5.1 General technical data

Length x Height x Width (without fiber) Laser Module Controller CKL 10xy z Controller CKL 24xy z SwitchBox SSL 1000	102.3 mm x 39 mm x 39 mm 145 mm x 43 mm x 60 mm 164 mm x 45 mm x 130 mm 67.2 mm x 66.2 mm x 28 mm
Weight Laser Module Controller CKL 10xy z Controller CKL 24xy z SwitchBox SSL 1000	approx. 0.7 kg approx. 0.29 kg approx. 0.35 kg approx. 0.2 kg
Operating voltage	100 ... 240 VAC / 50 ... 60 Hz
Power consumption	< 70 W
Maximum starting current	2.5 A
Operating environment	To be used in closed rooms only.
Protection Class	I
Degree of protection	IP 20
Overvoltage category	2
RFI suppression	DIN EN 61000-6-3:2007-09
Noise immunity	DIN EN 61000-6-2:2006-03
Ambient conditions Temperature range: Operation: Storage: Relative humidity: Operation: Storage: Altitude: Operation: Storage:	10 ... 40 °C -10 ... 60 °C ≤ 80 % ≤ 99 % 2,000 m 12,000 m
Mounting position	any

5.2 Laser types and useable accessories

Model name *LK www Tyzz

The indication: *, w, x, y, z are wildcards for following combinations:

* is a capital letter indicating the wavelength range of the laser	Letter	Wavelength [nm]
	U	300-399
	V	400-450
	B	451-500
	G	500-550
	Y	551-599
	R	600-699
	S	700-799
	T	800-899
	W	900-999
	Z	1000-1099
ww can be 00 to 99	indicating last two digits of the wavelength in nm	
x can be 0 to 499	indicating the optical output in mW	
y can be blank or S	indicates presence of a mechanical shutter (y=S)	
zz can be blank or any character	indicating client-specific, non safety relevant modifications	

Example

YLK 6150 TS DPSSL with wavelength 561nm and 50 mW output power and shutter

CKL 10xy z

OEM Controller with 12 V oder 24 V input voltage	
x=0 x>2	Standard model for 2 A drive current Drive current in A when anomalous from standard
y=0 y>2	Standard model for 2 V laser diode voltage Voltage in V when anomalous from standard
z	Digits, letters or series of digits and letters which show customized and not safety-related models

CKL 24xy z

Controller with 12 V oder 24 V input voltage	
x=0 x>2	Standard model for 2 A drive current Drive current in A when anomalous from standard
y=0 y>2	Standard model for 2 V laser diode voltage Voltage in V when anomalous from standard
z	Digits, letters or series of digits and letters which show customized and not safety-related models

Corresponding components:

DPSSL	Approved accessories	
*LK www Tyzz	Controller CKL 1000 SwitchBox SSL 1000 Power supply AEB70US24 Connection cable CKL1000 Power cord (purchased parts package of LASOS)	577009-4313-000 577009-4317-000 577001-9000-839 577007-4313-003
	Controller CKL 2400 or controller CKL 2400 OEM Power supply AEB70US24 Connection cable CKL 2400 Power cord (purchased parts package of LASOS)	577009-4307-000 577009-4306-000 577001-9000-839

Stand alone design generally consists of following components:

- DPSSL: *LK www TSzz (with shutter)
- Controller CKL 10xy z + SwitchBox SSL 1000
- Power supply (purchased parts package of LASOS)
- Cable (purchased parts package of LASOS)

or


- DPSSL: *LK www TSzz (with shutter)
- Controller CKL 24xy z
- Power supply (purchased parts package of LASOS)
- Cable (purchased parts package of LASOS)

DPSSL without shutter (*LK www Tzz) alone or in combination with Controller and/or SwitchBox and/or power supply are OEM versions. In this case the user is responsible for compliance with the corresponding safety regulations.


5.3 Dimension drawings:

DPSSL




STANDARD

OEM


SwitchBox SSL 1000


Controller


CKL 10xy z

CKL 24xy z

5.4 Maximum thermal output of the base plate of the laser module

While in operation, the temperature of the base plate of the laser module must not exceed 45 °C. Insufficient heat dissipation may cause overheating and total failure of the laser module.

6 LASOS Communicator

6.1 General

The LASOS Communicator allows operating laser models of the LASOS DPSSL series and LASOS LDM series. Basic information on model number and manufacturing information can be retrieved. The software gives information on the operation state of the laser and allows tuning the optical output power as well as switching On and Off.

Not all functions are supported by all laser models.

For the LASOS LDM series the software is only usable when the laser module is ordered with the stand alone option. The ComBox CDL 1000 is required. For details refer to the manual of the laser.

Product	USB	RS232
LASOS LDM series with ComBox CLD 1000 Models *LD xxxx and *LD xxxx F	Connector at ComBox CLD 1000	Connector at Multi-IO interface of ComBox CLD 1000
LASOS DPSS laser series with controller CKL 1000 Models *LK xxxx	Available at controller CKL 1000	Available with USB-RS232 adapter cable (optional)
LASOS DPSS laser series with controller CKL 2400 Models *LK xxxx	Not available	Available at controller CKL 2400

Table 1 Overview on availability of interfaces for the various laser models

6.2 Getting started

6.2.1 Installation of the software

The software is available as a Windows installer package in the folder "LASOS communicator" on the installation CD that is delivered with the laser. Double click on the file "LasosCommunicatorSetup.msi" starts the installation procedure automatically. The installation routine generates a folder "LASOS" in the Start Menu and an icon on the desktop.

Fig. 23 Desktop icon

6.2.2 Connection to the laser system

For availability of the interfaces refer to Table 1 of this manual and the manual of the particular laser model

6.2.2.1 Connecting by USB

Connect the laser system and the USB interface by a standard USB cable. When the laser module is connected for the first time the USB- driver has to be installed. Windows will prompt you to specify the location of the driver. Follow the instructions. The driver can be found in the folder "USB-drivers" on the installation CD.

When laser and computer are connected by USB a new serial port is created. This port has to be chosen to communicate with the laser (see chapter "Using the LASOS Communicator")

These instructions are also valid when an USB-RS232 adapter is used.

6.2.2.2 Connecting by RS232

Connect the RS232 connector of the computer with the interface connector of the laser system. Refer to the manual of the laser for the location of the interface (see also Table 1).

The following basic settings have to be applied.

Baud rate: 19,200

Data bit: 8

Parity: none

Stop bit: 1

Hand shake: none

6.3 Using the LASOS Communicator

Start the program by double click on the icon. The graphic user interface appears.






Fig. 24 Graphic user interface of the LASOS communicator

<p>Serial Communication</p> <p>Com Port: <input type="text" value="COM9"/></p> <p>Disconnect</p>	<p>If the laser is connected by RS232 choose the appropriate port number in the box.</p> <p>If the laser is connected by USB choose the appropriate port number that has been created after installing the drivers.</p> <p>The number can be chosen between 1 and 16. A free RS 232 port has to be available in this range.</p> <p>When the numbers are correct press the Connect/Disconnect button.</p>
--	--

Table 2 Functions of the LASOS communicator

	<p>The information fields specify the following parameters:</p> <ul style="list-style-type: none"> Diode current Output power Diode temperature Resonator temperature (DPSSL series only) Noise (DPSSL series only) Operation hours (DPSSL series only)
	<p>The middle upper part of the screen specifies the model number. The button Laser On/Off switches the laser on and off. The indicator "Warmup" blinks when the laser is in the warm up phase. The indicator "Error" lights when an error has occurred.</p>
	<p>The output power can be tuned by either by the slider or by typing a power value in the field below the slider. You can choose an output power between 0 and the nominal value (refer to the laser data sheet).</p>
	<p>In the right part of the screen additional information on the laser is shown:</p> <ul style="list-style-type: none"> Serial number (for reference) Shipping date (for reference) Wavelength Interlock open/closed Key switch on/off Laser status (On or Off)
	<p>The lower part of the screen shows a graph with the temporal behavior of the output power.</p>
	<p>With the button "Status report" you can save a file with status data of the laser. This can be sent to the manufacturer by e-mail for analysis in case of a failure. The information is encoded in the file and can only be retrieved by the manufacturer.</p> <p>The "Exit" button ends the program.</p>

When the diode temperature exceeds 40°C the error indicator lights and the laser is switched off. Check if the heat sink has sufficient capacity. The laser can be switched on again only after the electrical power has completely been disconnected and connected again.

EC DECLARATION OF CONFORMITY

according to Electromagnetic Compatibility Directive 2004/108/EC
of 2004/12/15

Herewith we declare that the following described product in his conception and design as well as in its construction we put into circulation is conform to the essential safety and healthiness requirements of the EMC Directive. In case of any modification not authorized by us regarding the product or the combination of laser-SwitchBox-Controller-power supply (incl. cables) this declaration loses its validity.

Manufacturer: LASOS Lasertechnik GmbH
Franz-Loewen-Str. 2
07745 Jena
Germany
Tel.: +49(0)3641 / 2944-0
www.lasos.com
info@lasos.com

Description of the product:

- Product name: LASOS® DPSS laser series
- Type: *LK www Tyzz (see the attachment of the EC declaration of conformity)
- Conventional use: for stand-alone and OEM-applications (e.g. flow cytometry, digital imaging, microscopy, measuring systems)

Technical documentation:

- Lists of parts, assembling plans, circuit diagrams, switching parts lists in the processing folder
- Example YLK 6120 T: Part list 577009-4186-000, date 2012-05-14
Assembling drawing 600009-4136-000, version 0A
Assembling drawing 600008-4136-000, version 00
Data sheet 577099-4186-000, version 0E
Data sheet 577099-4307-000, version 0B
Data sheet 577099-4313-000, version 0C
- Documentation of risk assessment
- Report on investigation of standards
- EMC-Test report

We declare the compliance with following standards which apply for the product:

- EN 61010-1:2011
Safety requirements for electrical equipment for measurement, control and laboratory use -
Part 1: General requirements

Used harmonized standards:

- EN 61000-6-1:2007
Electromagnetic compatibility (EMC) - Part 6-1: Generic standards - Immunity for residential, commercial and light-industrial environments
- EN 61000-6-2:2006
Electromagnetic compatibility (EMC) - Part 6-2: Generic standards - Immunity for industrial environments
- EN 61000-6-3:2011
Electromagnetic compatibility (EMC) - Part 6-3: Generic standards - Emission standard for residential, commercial and light-industrial environments
- EN 61000-6-4:2011
Electromagnetic compatibility (EMC) - Part 6-4: Generic standards - Emission standard for industrial environments

CE-documentation agent

Name: Rainer Laske, Designing Engineer

Information/Identity about the signatory:

Dr. Meyer, Managing Director

Date / Signature:.....

11.2.13

Attachment of EC Declaration of Conformity

Model name *LK www Tyzz

The indication: *, w, x, y, z are wildcards for following combinations:

* is a capital letter indicating the wavelength range of the laser	Letter	Wavelength [nm]
	U	300-399
	V	400-450
	B	451-500
	G	500-550
	Y	551-599
	R	600-699
	S	700-799
	T	800-899
	W	900-999
	Z	1000-1099
ww can be 00 to 99	indicating last two digits of the wavelength in nm	
x can be 0 to 499	indicating the optical output in mW	
y can be blank or S	indicates presence of a mechanical shutter (y=S)	
zz can be blank or any character	indicating client-specific, non safety relevant modifications	

Example

YLK 6150 TS DPSSL with wavelength 561nm and 50 mW output power and shutter

CKL 10xy z

OEM Controller with 12 V oder 24 V input voltage	
x=0 x>2	Standard model for 2 A drive current Drive current in A when anomalous from standard
y=0 y>2	Standard model for 2 V laser diode voltage Voltage in V when anomalous from standard
z	Digits, letters or series of digits and letters which show customized and not safety-related models

CKL 24xy z

Controller with 12 V oder 24 V input voltage	
x=0 x>2	Standard model for 2 A drive current Drive current in A when anomalous from standard
y=0 y>2	Standard model for 2 V laser diode voltage Voltage in V when anomalous from standard
z	Digits, letters or series of digits and letters which show customized and not safety-related models

Corresponding components:

DPSSL	Approved accessories	
*LK www Tyzz	Controller CKL 1000 SwitchBox SSL 1000 Power supply AEB70US24 Connection cable CKL1000 Power cord (purchased parts package of LASOS)	577009-4313-000 577009-4317-000 577001-9000-839 577007-4313-003
	Controller CKL 2400 or controller CKL 2400 OEM Power supply AEB70US24 Connection cable CKL 2400 Power cord (purchased parts package of LASOS)	577009-4307-000 577009-4306-000 577001-9000-839

Stand alone design generally consists of following components:

- DPSSL: *LK www TSzz (with shutter)
- Controller CKL 10xy z + SwitchBox SSL 1000
- Power supply (purchased parts package of LASOS)
- Cable (purchased parts package of LASOS)

or

- DPSSL: *LK www TSzz (with shutter)
- Controller CKL 24xy z
- Power supply (purchased parts package of LASOS)
- Cable (purchased parts package of LASOS)

DPSSL without shutter (*LK www Tzz) alone or in combination with Controller and/or SwitchBox and/or power supply are OEM versions. In this case the user is responsible for compliance with the corresponding safety regulations.

EC DECLARATION OF CONFORMITY

according to Low Voltage Directive 2006/95/EC
of 2006/12/12

Herewith we declare that the following described product in his conception and design as well as in its construction we put into circulation is conform to the essential safety and healthiness requirements of the Low Voltage Directive. In case of any modification not authorized by us regarding the product or the combination of laser-SwitchBox-Controller-power supply (incl. cables) this declaration loses its validity.

Manufacturer: LASOS Lasertechnik GmbH
Franz-Loewen-Str. 2
07745 Jena
Germany
Tel.: +49(0)3641 / 2944-0
www.lasos.com
info@lasos.com

Description of the product:

- Product name: LASOS® DPSS laser series
- Types: *LK wwx TSzz (see the attachment of the EC declaration of conformity)
- Conventional use: for stand-alone and OEM-applications (e.g. flow cytometry, digital imaging, microscopy, measuring systems)

Technical documentation:

- Lists of parts, assembling plans, circuit diagrams, switching parts lists in the processing folder
- Example YLK 6120 T: Part list 577009-4186-000, date 2012-05-14
Assembling drawing 600009-4136-000, version 0A
Assembling drawing 600008-4136-000, version 00
Data sheet 577099-4186-000, version 0E
Data sheet 577099-4307-000, version 0B
Data sheet 577099-4313-000, version 0C
- Documentation of risk assessment
- Report on investigation of standards
- EMC-Test report

It will declare the compliance with further directives which are applying for the product:

- EMC-directive 2004/108/EC

Used harmonized standards:

- EN 60825-1:2008
Safety of laser products - Part 1: Equipment classification and requirements
- EN 61000-6-1:2007
Electromagnetic compatibility (EMC) - Part 6-1: Generic standards - Immunity for residential, commercial and light-industrial environments
- EN 61000-6-2:2006
Electromagnetic compatibility (EMC) - Part 6-2: Generic standards - Immunity for industrial environments
- EN 61000-6-3:2011
Electromagnetic compatibility (EMC) - Part 6-3: Generic standards - Emission standard for residential, commercial and light-industrial environments
- EN 61000-6-4:2011
Electromagnetic compatibility (EMC) - Part 6-4: Generic standards - Emission standard for industrial environments
- EN 61010-1:2011
Safety requirements for electrical equipment for measurement, control and laboratory use - Part 1: General requirements

CE-documentation agent

Name: Rainer Laske, Designing Engineer

Information/Identity of the signatory:

Dr. Meyer, Managing Director

Date / Signature:

19.2.13

Attachment of EC Declaration of Conformity

Model name ***LK ww x TSzz**

The indication: *, w, x, y, z are wildcards for following combinations:

* is a capital letter indicating the wavelength range of the laser	Letter	Wavelength [nm]
	U	300-399
	V	400-450
	B	451-500
	G	500-550
	Y	551-599
	R	600-699
	S	700-799
	T	800-899
	W	900-999
	Z	1000-1099
ww can be 00 to 99	indicating last two digits of the wavelength in nm	
x can be 0 to 499	indicating the optical output in mW	
y can be blank or S	indicates presence of a mechanical shutter (y=S)	
zz can be blank or any character	indicating client-specific, non safety relevant modifications	

Example

YLK 6150 TS

DPSSL with wavelength 561nm and 50 mW output power and shutter

CKL 10xy z

OEM Controller with 12 V oder 24 V input voltage	
x=0 x>2	Standard model for 2 A drive current Drive current in A when anomalous from standard
y=0 y>2	Standard model for 2 V laser diode voltage Voltage in V when anomalous from standard
z	Digits, letters or series of digits and letters which show customized and not safety-related models

CKL 24xy z

Controller with 12 V oder 24 V input voltage	
x=0 x>2	Standard model for 2 A drive current Drive current in A when anomalous from standard
y=0 y>2	Standard model for 2 V laser diode voltage Voltage in V when anomalous from standard
z	Digits, letters or series of digits and letters which show customized and not safety-related models

Corresponding components:

DPSSL	Approved accessories	
*LK wwx TSzz	Controller CKL 1000 SwitchBox SSL 1000 Power supply AEB70US24 Connection cable CKL1000 Power cord (purchased parts package of LASOS)	577009-4313-000 577009-4317-000 577001-9000-839 577007-4313-003
	Controller CKL 2400 or controller CKL 2400 OEM Power supply AEB70US24 Connection cable CKL 2400 Power cord (purchased parts package of LASOS)	577009-4307-000 577009-4306-000 577001-9000-839

Stand alone design generally consists of following components:

- DPSSL: *LK wwx TSzz (with shutter)
- Controller CKL 10xy z + SwitchBox SSL 1000
- Power supply (purchased parts package of LASOS)
- Cable (purchased parts package of LASOS)

or

- DPSSL: *LK wwx TSzz (with shutter)
- Controller CKL 24xy z
- Power supply (purchased parts package of LASOS)
- Cable (purchased parts package of LASOS)

DPSSL without shutter (*LK wwx Tzz) alone or in combination with Controller and/or SwitchBox and/or power supply are OEM versions. In this case the user is responsible for compliance with the corresponding safety regulations.